



## PRODUCT SPECIFICATION

catalog number:	<b>A60123</b>
product description:	Influenza B Monoclonal Antibody
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	20 mM PBS, 50 mM NaCl, pH 7.4
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG <sub>2a</sub>
concentration:	[lot specific] mg/mL
specificity:	Specific for Influenza B virus. Recognizes the nucleoprotein of numerous strains.
cross-reactivity:	No cross-reactivity with Parainfluenza, RSV, Adenovirus.
pairing:	This antibody can be used as labeling with A60124574P.
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.