



PRODUCT SPECIFICATION

catalog number:	A60130
product description:	Influenza B Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01M PBS, pH 7.2
preservative:	0.1% NaN ₃
isotype:	IgG _{2b}
concentration:	[lot specific] mg/mL
specificity:	Specific for nuclear protein
pairing:	This antibody pairs well with A60110 (used as the conjugate).
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.