



## PRODUCT SPECIFICATION

catalog number:	<b>A60223</b>
product description:	Parainfluenza Virus Type 3 Monoclonal Antibody
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01 M PBS, pH 7.2
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG2a
concentration:	[lot specific] mg/mL
specificity:	Specific for hemagglutinin of Parainfluenza virus. Positive for Bovine PI3.
pairing:	This antibody can be paired with A60224 in ELISA or IFA.
notes/applications:	ELISA and IFA
expiration:	2 years from date of receipt when stored at 2-8°C
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.