



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | A60225 |
| product description: | Parainfluenza Virus Type 1 Monoclonal Antibody |
| lot number: | [lot specific] |
| source: | Murine ascites |
| form: | Liquid |
| purification method: | Protein A chromatography |
| purity: | >90% |
| buffer: | 0.01 M PBS, pH 7.2 |
| preservative: | 0.1% NaN ₃ |
| isotype: | IgG2b |
| concentration: | [lot specific] mg/mL |
| specificity: | Specific for Parainfluenza virus, type 1. Negative against types 2 & 3 by indirect IFA. Reactive with Sendai virus. |
| notes/applications: | ELISA and IFA |
| expiration: | 2 years from date of receipt when stored at 2-8°C |
| storage: | 2-8°C |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |