



PRODUCT SPECIFICATION

catalog number:	A60226
product description:	Parainfluenza Virus Type 1 Monoclonal Antibody
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01 M PBS, pH 7.2
preservative:	0.1% NaN ₃
isotype:	IgG2b
concentration:	[lot specific] mg/mL
specificity:	Specific for Parainfluenza virus, type 1. Negative against types 2 & 3 by indirect IFA. Reactive with Sendai virus.
notes/applications:	ELISA and IFA
expiration:	2 years from date of receipt when stored at 2-8°C
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.