



## PRODUCT SPECIFICATION

catalog number:	<b>A60234</b>
product description:	<i>Campylobacter jejuni</i> Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Total <i>C. jejuni</i> sonicate
source:	Murine; in vitro production
form:	Liquid
purification method:	Protein A/G chromatography
buffer:	10mM PBS, pH 7.4
preservative:	0.05% NaN <sub>3</sub>
isotype:	IgM
concentration:	[lot specific] mg/mL
specificity:	This monoclonal is specific for the non-flagellar antigen of both <i>C. jejuni</i> types 1 and 2.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.