



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A60250**

PRODUCT:	Parainfluenza III Monoclonal Antibody
FORM:	Purified Immunoglobulin
BUFFER:	0.02M PBS, 0.25M NaCl, pH 7.6
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG <sub>1</sub>
STORAGE CONDITIONS:	2-8° C
SPECIFICITY:	Specific for the HN protein common to all Parainfluenza III strains. Exhibits no cross-reactivity to Parainfluenza I or II.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.