



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A60332**

PRODUCT:	Mycoplasma pneumoniae Monoclonal Antibody
LOT NUMBER:	[lot specific]
FORM:	Immunoglobulin purified from mouse ascites fluid or culture medium
PURIFICATION:	Protein A chromatography
PURITY:	>90%
BUFFER:	0.01 M phosphate buffered saline, pH 7.2
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG <sub>1</sub>
CONCENTRATION:	[lot specific]
STORAGE CONDITIONS:	Short term: 2-8°C; Long term: -20 °C
SPECIFICITY:	Reactive with a >200 kDa antigen of Mycoplasma pneumoniae.
APPLICATION:	ELISA, IFA & WB. Neg. with uninfected cells
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.