



## PRODUCT SPECIFICATION

catalog number:	<b>A60342</b>
product description:	Mycoplasma Pneumoniae P1
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01 M PBS, pH 7.2
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG2b
concentration:	[lot specific] mg/mL
specificity:	Specific for the P1 adhesin of M. pneumoniae
pairing:	Use as capture; pair with A60345 as conjugate.
notes/applications:	ELISA and IFA
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.