



PRODUCT SPECIFICATION

catalog number:	A60401
product description:	<i>Listeria monocytogenes</i> Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>90% pure
buffer:	0.01M PBS, pH 7.2
preservative:	0.1% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	This monoclonal antibody reacts with <i>L. monocytogenes</i> .
pairing:	This antibody can be used as the capture antibody along with A60402041P.
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.