



## PRODUCT SPECIFICATION

catalog number:	<b>A60402</b>
product description:	<i>Listeria monocytogenes</i> Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>90% pure
buffer:	0.01M PBS, pH 7.2
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG <sub>2a</sub>
concentration:	[lot specific] mg/mL
specificity:	This monoclonal antibody reacts with a membrane antigen of <i>L. monocytogenes</i> .
pairing:	This antibody can be used as the conjugate with A60401041P.
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.