



PRODUCT SPECIFICATION

catalog number:	A60403
product description:	<i>Listeria monocytogenes</i> Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>90% pure
buffer:	0.01M PBS, pH 7.2
preservative:	0.1% NaN ₃
isotype:	IgG _{2b}
concentration:	[lot specific] mg/mL
specificity:	This monoclonal antibody is specific for the p60 secreted protein of <i>L. monocytogenes</i> .
pairing:	This antibody can be used as the capture with A60404041P or A60406041P; this can also be used as the conjugate with A60404401P or A60405401P
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.