



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A60520**

PRODUCT:	<u>Treponema pallidum</u> Monoclonal Antibody
IMMUNOGEN:	<u>Treponema pallidum</u>
PURIFICATION:	Immunoglobulin purified on protein A sepharose
BUFFER:	PBS pH 7.4
PRESERVATIVE:	0.1% NaN <sub>3</sub>
SUBCLASS:	IgG1
STORAGE CONDITIONS:	2°C - 8°C
REACTIVITY:	Reacts with ultrasonicated lysates of <u>Tr. pallidum</u> and <u>Tr. reiter</u> .
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.