



PRODUCT SPECIFICATION

catalog number:	A60530
product description:	<i>Salmonella typhimurium</i> Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Lipopolysaccharides (LPS) of <i>S. typhimurium</i>
source:	Murine
form:	Purified liquid
purification method:	Protein A chromatography
buffer:	PBS, pH 7.4
preservative:	0.1% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	LPS of <i>S. typhimurium</i>
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.