



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | A60575 |
| product description: | <i>Salmonella</i> (A Group) Monoclonal Antibody, Purified |
| lot number: | [lot specific] |
| immunogen: | Heat inactivated cells from <i>S. paratyphi</i> , <i>S. typhimurium</i> and <i>S. enteritidis</i> |
| source: | Murine |
| form: | Liquid |
| purification method: | Protein G chromatography |
| buffer: | PBS, pH 7.4 |
| preservative: | 0.1% NaN ₃ |
| isotype: | IgG ₁ |
| concentration: | [lot specific] mg/mL |
| specificity: | Specific for sero group A and reacting with the 0-2 Ag determinant, this monoclonal is known to cross-react with <i>S. paratyphi</i> . |
| expiration date: | [lot specific] |
| storage: | 2-8°C |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |