



PRODUCT SPECIFICATION

catalog number:	A60580
product description:	Pan <i>Salmonella</i> Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Heat inactivated cells from <i>S. paratyphi</i> , <i>S. typhimurium</i> and <i>S. enteritidis</i>
source:	Murine
form:	Liquid
purification method:	Protein G chromatography
buffer:	PBS, pH 7.4
preservative:	0.1% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	This antibody is pan specific to <i>Salmonella reacting with the core Ag</i> determinant. It is known to cross-react with <i>S. paratyphi</i> , <i>S. typhimurium</i> , <i>S. enteritidis</i> , <i>S. choleraesuis</i> , <i>S. newport</i> , <i>S. anatum</i> and <i>S. selandia</i> .
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.