



PRODUCT SPECIFICATION

CATALOG NUMBER: A60600

PRODUCT:	<u>Escherichia coli</u> Monoclonal Antibody
IMMUNOGEN:	Purified <u>E. coli</u> serotype 0157
PURITY:	Purified by Protein A
BUFFER:	10 mM Phosphate, pH 7.4, 150 mM NaCl
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG3
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	<u>E. coli</u> LPS and O polysaccharide