



## PRODUCT SPECIFICATION

catalog number:	<b>A60650</b>
product description:	MRSA Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)
source:	Murine; in vitro production
form:	Liquid
purification method:	Protein A/G chromatography
buffer:	10mM PBS, pH 7.4
preservative:	0.05% NaN <sub>3</sub>
isotype:	IgG <sub>2a</sub> , kappa
concentration:	[lot specific] mg/mL
specificity:	This antibody is specific for MRSA also known as ORSA (oxacillin-resistant <i>Staph aureus</i> ).
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.