



## PRODUCT SPECIFICATION

|                      |   |
|----------------------|---|
| catalog number:      | <b>A60695</b>   |
| product description: | Norovirus Genogroup I (GI) Monoclonal Antibody, Purified  |
| lot number:          | [lot specific]  |
| source:              | Murine ascites  |
| form:                | Liquid  |
| purification method: | Protein A chromatography  |
| purity:              | > 90% pure  |
| buffer:              | 0.01M PBS, pH 7.2   |
| preservative:        | 0.1% NaN <sub>3</sub>   |
| isotype:             | IgG <sub>2a</sub>   |
| concentration:       | [lot specific] mg/mL  |
| specificity:         | This monoclonal antibody recognizes the GI capsid protein.  |
| cross-reactivity:    | It does not cross-react with the GII capsid protein.  |
| pairing:             | This antibody can be used as the capture antibody and pairs with A60691 as the detection antibody.  |
| storage:             | 2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.   |
| hazard/biohazard:    | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |