



## PRODUCT SPECIFICATION

catalog number:	<b>A60720</b>
product description:	Norovirus Genogroup II (GII) Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	> 90% pure
buffer:	0.01M PBS, pH 7.2
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
specificity:	This monoclonal antibody recognizes the GI capsid protein.
cross-reactivity:	It does not cross-react with the GI capsid protein.
pairing:	This antibody can be used as the conjugate antibody in a sandwich EIA and pairs with A60740 as the capture antibody.
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.