



PRODUCT SPECIFICATION

catalog number:	A60750
product description:	Smallpox B5R Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Recombinant B5R protein from <i>Variola</i> virus
source:	Murine; in vitro production
form:	Liquid
purification method:	Protein A/G chromatography
buffer:	10mM PBS, pH 7.4
preservative:	0.05% NaN ₃
isotype:	IgG ₁ , kappa
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.