



PRODUCT SPECIFICATION

catalog number:	A60830
product description:	<i>Pseudomonas aeruginosa</i> (Serotype 6C) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	<i>P. aeruginosa</i> serotype 6C
source:	Murine
form:	Liquid
purification method:	Protein A/G chromatography
buffer:	10mM PBS, pH 7.4
preservative:	0.05% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	This antibody is specific for serotype 6C of <i>S. aeruginosa</i> .
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.