



## PRODUCT SPECIFICATION

catalog number:	<b>A61143</b>
product description:	Insulin Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Recombinant insulin
source:	Murine ascites
form:	Colorless liquid
purification method:	Affinity chromatography
buffer:	15mM phosphate buffer, 0.15M NaCl, pH 7.2
preservative:	0.09% NaN <sub>3</sub>
isotype:	IgG <sub>1k</sub>
concentration:	[lot specific] mg/mL ( $A^{1\%}_{280nm} = 14.0$ )
specificity:	Human insulin
application/notes:	Can be used in a sandwich ELISA as the detection antibody in combination with A61142 used as the capture antibody.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.