



PRODUCT SPECIFICATION

catalog number:	A61169
product description:	FITC Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	FITC-BSA
source:	Murine
form:	Liquid
purity:	>95% by SDS-PAGE
buffer:	PBS, pH 7.0
preservative:	0.02% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL (based on A280 _{nm} , 1 mg/ml, 1cm= 1.4)
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
date of manufacture:	[lot specific]
expiry:	12 months from date of receipt when stored at -20°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.