



PRODUCT SPECIFICATION

CATALOG NUMBER: A61210

PRODUCT:	IGF-1 Monoclonal Antibody
SOURCE:	Tissue culture supernatant
IMMUNOGEN:	IGF-1-Ovalbumin
PURIFICATION METHOD:	Protein A affinity chromatography
BUFFER:	PBS containing 1 mg/mL BSA
PRESERVATIVE:	10mM NaN ₃
SUBCLASS:	IgG ₁
AFFINITY CONSTANT:	1 x 10 ¹⁰ L/mole
STORAGE CONDITIONS:	2-8°C short term, -20°C long term
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.