



PRODUCT SPECIFICATION

CATALOG NUMBER: A61230

PRODUCT:	IGF-2 Monoclonal Antibody
SOURCE:	Mouse ascites fluid
PURIFICATION METHOD:	Ion exchange chromatography
BUFFER:	PBS
PRESERVATIVE:	10mM NaN ₃
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	4-8°C short term, -20°C long term
CROSS REACTIVITY:	IGF-1 <0.01% and insulin <0.01%
APPLICATIONS:	Measurement of serum IGF-2 levels. Application may be RAI, ELISA, IRMA or immunohistochemistry.
HAZARD/BIOHAZARD:	This material contains sodium azide as preservative, take appropriate precaution. There is no known biohazard associated with this product.