



PRODUCT SPECIFICATION

CATALOG NUMBER: A61290

PRODUCT:	Human Insulin Monoclonal Antibody
SOURCE:	Immunoglobulin from mouse ascites
IMMUNOGEN:	Human Insulin
PURIFICATION:	Purified by chromatography on protein A Sepharose. Purity is tested by electrophoresis.
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG ₁ , kappa
STORAGE CONDITIONS:	4°C
REACTIVITY:	These antibodies cross-react with human proinsulin, bovine insulin (30%) and porcine insulin.
APPLICATION:	Detection of insulin in two-site enzyme immunoassay.