



## PRODUCT SPECIFICATION

catalog number:	<b>A61310087P</b>
product description:	C-Peptide Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human c-peptide
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
buffer:	PBS, pH 7.2
preservative:	0.05% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
application/notes:	RIA, EIA, immunohistochemistry
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.