



PRODUCT SPECIFICATION

catalog number:	A61370087P
product description:	C-Peptide Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human c-peptide
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
purity:	>90% by SDS-PAGE
buffer:	PBS, pH 7.2
preservative:	0.05% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
application/notes:	RIA, EIA, immunohistochemistry
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.