



## PRODUCT SPECIFICATION

catalog number:	<b>A61372</b>
product description:	C-Peptide Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	C-Peptide conjugate BSA
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
purity:	>98%
buffer:	50 mM NaPO <sub>4</sub> Buffer with 0.9% NaCl, pH 7.2±0.2, 0.02% Proclin
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
application/notes:	RIA, EIA, immunohistochemistry, and western blot
storage:	2-8°C
hazard/biohazard:	There is no known biohazard associated with this product.