



## PRODUCT SPECIFICATION

catalog number:	<b>A61405</b>
product description:	IGFBP-3 Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Human Recombinant IGFBP-3
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	BBS 0.1M pH 8.0
preservative:	none
isotype:	N/A
concentration:	[lot specific] mg/mL
specificity:	Recognizes natural and recombinant human IGFBP-3
pairing:	This antibody can be used in sandwich ELISA as the conjugate antibody in conjugation with A61406 as the capture
storage:	-20°C
hazard/biohazard:	There is no known biohazard associated with this product.