



PRODUCT SPECIFICATION

CATALOG NUMBER: A61431

PRODUCT:	Glucagon Monoclonal Antibody
IMMUNOGEN:	Human recombinant Glucagon
FORM:	Liquid
PURIFICATION METHOD:	Protein A affinity chromatography
BUFFER:	PBS, pH 7.2
PRESERVATIVE:	0.05% sodium azide
SUBCLASS:	IgG ₁
BINDING CONSTANT:	5×10^{10}
STORAGE CONDITIONS:	2-8°C
APPLICATIONS:	RIA, EIA, Immunochemistry
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.