



PRODUCT SPECIFICATION

catalog number:	A62014
product description:	Adenovirus Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	> 90% pure
buffer:	0.01M PBS, pH 7.2
preservative:	0.1% NaN ₃
isotype:	IgG _{2a}
concentration:	[lot specific] mg/mL
specificity:	This monoclonal antibody is specific for the hexon group antigen of numerous adenovirus; known reactivity with at least 16 serotypes including types 40 and 41.
pairing:	This antibody pairs well in a sandwich EIA with A62010.
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.