



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A62050**

PRODUCT:	Adenovirus Monoclonal Antibody
IMMUNOGEN:	Adeno 3
FORM:	Purified immunoglobulin from ascites
BUFFER:	0.02M phosphate buffer, 0.25M sodium chloride, pH 7.6
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG <sub>1k</sub>
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	All 41 serotypes of Adenovirus
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.