



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: A62060**

PRODUCT:	Adenovirus Hexon Monoclonal Antibody
SOURCE:	Immunoglobulin from mouse ascites
IMMUNOGEN:	Hexon antigen of human adenovirus
PURIFICATION:	Purified by chromatography on protein A Sepharose. Purity is tested by electrophoresis.
FORM:	PBS, pH 7.4
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG <sub>2a</sub>
STORAGE CONDITIONS:	4°C
APPLICATION:	Can be used in ELISA and immunodiffusion.