



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A62070**

PRODUCT:	Adenovirus Monoclonal Antibody
PURIFICATION:	Purified by chromatography on Protein A Sepharose. Purity is tested by electrophoresis.
FORM:	Immunoglobulin
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% of sodium azide
SUBCLASS:	IgG <sub>2a</sub>
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Hexon antigen of human adenovirus
APPLICATION:	ELISA