



PRODUCT SPECIFICATION

CATALOG NUMBER: A63221

PRODUCT:	Clostridium botulinum D toxoid Polyclonal
SOURCE:	Ascites fluid
PURIFICATION:	Protein A chromatography
PURITY:	>90%
BUFFER:	0.01 M PBS, pH 7.2
PRESERVATIVE:	0.1% sodium azide
STORAGE CONDITIONS:	-20°C until ready to use. Avoid repeated freeze-thawing cycles.
SPECIFICITY:	Specific for toxin/toxoid D of C. botulinum.
APPLICATION:	Functions in ELISA and IFA.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.