



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A63230**

PRODUCT:	Clostridium Difficile Toxin A Monoclonal Antibody
SOURCE:	Ascites fluid
PURIFICATION:	Purified from ascites fluid by Protein A chromatography
PURITY:	>90%
BUFFER:	0.01 M PBS, pH 7.2
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG <sub>1</sub>
STORAGE CONDITIONS:	-20°C until ready to use. Avoid repeated freeze-thawing cycles.
SPECIFICITY:	Specific for the toxin/toxoid A of Clostridium difficile.
APPLICATION:	Negative versus other toxins. ELISA.
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.