



PRODUCT SPECIFICATION

catalog number:	A63280
product description:	<i>Clostridium difficile</i> Toxin B Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	> 90% pure
buffer:	0.01M PBS, pH 7.2
preservative:	0.1% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	This monoclonal antibody is specific for the toxin B of <i>Clostridium difficile</i> .
cross-reactivity:	There is minimal cross reaction with toxin A.
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.