



PRODUCT SPECIFICATION

CATALOG NUMBER: A63300

| | |
|----------------------|--|
| PRODUCT DESCRIPTION: | Clostridium difficile Toxin B Monoclonal Antibody |
| IMMUNOGEN: | <i>C.difficile</i> Toxin B |
| SOURCE: | Ascites fluid |
| PURIFICATION: | Protein A chromatography |
| PURITY: | >90% |
| BUFFER: | 0.01 M PBS, pH 7.2 |
| PRESERVATIVE: | 0.1% sodium azide |
| SPECIFICITY: | Specific for the toxin B of <i>Clostridium difficile</i> . Negative vs. other toxins |
| SUBCLASS: | IgG ₁ |
| STORAGE: | 2-8°C short term, -20°C long term Avoid repeated freeze-thaw cycles |
| HAZARD/BIOHAZARD: | This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product. |