



## PRODUCT SPECIFICATION

catalog number:	<b>A63360</b>
product description:	<i>Clostridium difficile</i> (Toxin A) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	<i>C. difficile</i> Toxin A
source:	Murine; in vitro production
form:	Liquid
purification method:	Protein A/G chromatography
buffer:	10mM PBS, pH 7.4
preservative:	0.05% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub> , kappa
molecular weight:	308kDa
concentration:	[lot specific] mg/mL
specificity:	This antibody is specific for <i>C. difficile</i> Toxin A
cross-reactivity:	This monoclonal does not react with Cholera subunit a, Cholera toxin, <i>P. aeruginosa</i> exotoxin A, H-LT or P-LT.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.