



PRODUCT SPECIFICATION

CATALOG NUMBER: A63530

PRODUCT:	Rotavirus Monoclonal Antibody
IMMUNOGEN:	Bovine Rotavirus isolate
FORM:	Purified immunoglobulin from mouse ascites
PURIFICATION:	Protein A Sepharose
BUFFER:	PBS, pH 7.2 with 0.01% sodium azide
SUBCLASS:	IgG _{2a}
STORAGE CONDITIONS:	2-8°C short term -20°C long term
APPLICATIONS:	ELISA and Western blot.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.