



PRODUCT SPECIFICATION

CATALOG NUMBER: A63550

PRODUCT DESCRIPTION:	Rotavirus Monoclonal Antibody
PURIFICATION:	Protein A chromatography
PURITY:	>90%
BUFFER:	0.01 M PBS, pH 7.2
PRESERVATIVE:	0.1% NaN ₃
SUBCLASS:	IgG _{2a}
SPECIFICITY:	Reactive with intact virus of strain, RRV, WA & bovine. Pairs with A63540 (capture antibody) in a lateral flow assay.
STORAGE:	-20°C until ready to use
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.