



## PRODUCT SPECIFICATION

catalog number:	<b>A63650</b>
product description:	Malaria <i>P. vivax</i> LDH Monoclonal Antibody, Purified
lot number:	[lot specific]
form:	Liquid
purification method:	Protein A chromatography
purity:	>95% by electrophoresis
buffer:	10mM PBS, pH7.4
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG <sub>2b</sub>
concentration:	[lot specific] mg/mL
specificity:	This monoclonal is recognizes pLDH from <i>P. vivax</i> .
pairing:	This antibody can be used as the capture antibody and paired with A63610 as the conjugate antibody.
expiration date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.