



## PRODUCT SPECIFICATION

catalog number:	<b>A63670</b>
product description:	Malaria P.v. pLDH Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	<i>Plasmodium vivax</i> LDH, recombinant
source:	Murine
form:	Liquid
purification method:	Protein G chromatography
purity:	>95%
buffer:	10mM PBS, 50mM NaCl, pH 8.0
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
specificity:	This monoclonal is specific to malaria <i>Plasmodium vivax</i> lactate dehydrogenase (Pv LDH)
pairing:	This antibody can be used in sandwich ELISA as the capture antibody in conjunction with A63660 (as the detection antibody).
expiration date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.