



## PRODUCT SPECIFICATION

catalog number:	<b>A63730</b>
product description:	Malaria HRP-II Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Recombinant malaria HRP-II
source:	Murine
form:	Liquid
purity:	>90% by SDS-PAGE and HPLC
buffer:	PBS, pH7.4
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
specificity:	This monoclonal is specific for the histidine rich protein-II (HRP-II) from <i>P. falciparum</i> .
pairing:	This antibody can be used as the conjugate or capture antibody and pairs with A63710 or A63720. It is known not to self-pair well.
expiration date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.