



PRODUCT SPECIFICATION

catalog number:	A63740
product description:	Malaria HRP-II Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Recombinant malaria HRP-II
source:	Murine
form:	Liquid
purity:	>95% by SDS-PAGE
buffer:	PBS, pH7.2
preservative:	0.02% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	This monoclonal is specific for <i>P. falciparum</i> histidine rich protein-II (HRP-II).
pairing:	This antibody can be used as the conjugate antibody and paired with A63750 as the capture antibody.
expiration date:	[lot specific]
storage:	-20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.