



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | A63830 |
| product description: | Malaria pf-HRP-2 Monoclonal Antibody, Purified |
| lot number: | [lot specific] |
| immunogen: | Native Plasmodium falciparum HRP-2 |
| source: | Murine |
| form: | Liquid |
| purity: | >90% by SDS-PAGE |
| buffer: | 0.15 Potassium Phosphate Buffer, 0.85% Sodium Chloride, pH 7.2 |
| preservative: | 0.05% NaN ₃ |
| isotype: | IgG ₁ , kappa |
| concentration: | [lot specific] mg/mL |
| specificity: | Plasmodium falciparum (malaria) HRP-2 |
| pairing: | This antibody forms a match pair with Malaria HRP-2 Recombinant Antigen J63710314. |
| notes/applications: | ELISA, Lateral Flow |
| storage: | 2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles. |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |