



PRODUCT SPECIFICATION

catalog number:	A63840
product description:	Malaria pf-HRP-2 Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Native Plasmodium falciparum HRP-2
source:	Murine
form:	Liquid
purity:	>90% by SDS-PAGE
buffer:	0.15M Potassium Phosphate Buffer, 0.85% Sodium Chloride, pH 7.2
preservative:	0.05% NaN ₃
isotype:	IgG _{1, kappa}
concentration:	[lot specific] mg/mL
specificity:	Plasmodium falciparum (malaria) HRP-2
pairing:	This antibody forms a match pair with Malaria HRP-2 Recombinant Antigen J63710314.
notes/applications:	ELISA, Lateral Flow
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.