



PRODUCT SPECIFICATION

catalog number:	A66012
product description:	ACTH (N-Terminal) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human ACTH N-terminal fragment
source:	Murine; in vitro production
form:	Liquid
purification method:	Protein A/G chromatography
buffer:	10mM PBS, pH 7.4
preservative:	0.05% NaN ₃
isotype:	IgG _{1, kappa}
concentration:	[lot specific] mg/mL
specificity:	This antibody is specific for human Synacthen (aa1-24 of ACTH). It does not react with CLIP (aa25-39 of ACTH).
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.